



ClearCanvas ImageServer

Dicom Conformance Statement

Document Version: 1.0

Product Name(s): ClearCanvas ImageServer

Release: 1.3

Date: December 11, 2008

1. CONFORMANCE STATEMENT OVERVIEW

The ClearCanvas ImageServer is an open source DICOM server. The ImageServer software is composed of a web based interface used for administration and a Windows service. The ImageServer has a scalable design that allows it to run on a single server or a cluster of servers.

The application supports long term storage of images, waveforms, reports, and measurements. It also supports querying of its contents by remote systems, and retrieval of the stored objects. It also contains a rules engine for defining retention times of studies, compression of online studies, and auto-routing of studies to other DICOM devices.

**Table 1-1
NETWORK SERVICES**

SOP Classes	User of Service (SCU)	Provider of Service (SCP)
Image Transfer		
Computed Radiography Image Storage	Yes	Yes
CT Image Storage	Yes	Yes
Digital Intra-oral X-Ray Image Storage – For Presentation	Yes	Yes
Digital Intra-oral X-Ray Image Storage – For Processing	Yes	Yes
Digital Mammography X-Ray Image Storage – For Presentation	Yes	Yes
Digital Mammography X-Ray Image Storage – For Processing	Yes	Yes
Digital X-Ray Image Storage – For Presentation	Yes	Yes
Digital X-Ray Image Storage – For Processing	Yes	Yes
Enhanced CT Image Storage	Yes	Yes
Enhanced MR Image Storage	Yes	Yes
Enhanced XA Image Storage	Yes	Yes
Enhanced XRF Image Storage	Yes	Yes
Hardcopy Grayscale Image Storage SOP Class (Retired)	Yes	Yes
Hardcopy Color Image Storage SOP Class (Retired)	Yes	Yes
MR Image Storage	Yes	Yes
Multi-frame 1 Color Secondary Capture Image Storage	Yes	Yes
Multi-frame Grayscale Byte Secondary Capture Image Storage	Yes	Yes
Multi-frame Grayscale Word Secondary Capture Image Storage	Yes	Yes
Multi-frame Single Bit Secondary Capture Image Storage	Yes	Yes
Nuclear Medicine Image Storage (Retired)	Yes	Yes
Nuclear Medicine Image Storage	Yes	Yes
Ophthalmic Photography 16 Bit Image Storage	Yes	Yes
Ophthalmic Photography 8 Bit Image Storage	Yes	Yes
Ophthalmic Tomography Image Storage	Yes	Yes
Positron Emission Tomography Image Storage	Yes	Yes

RT Image Storage	Yes	Yes
Secondary Capture Image Storage	Yes	Yes
Ultrasound Image Storage	Yes	Yes
Ultrasound Image Storage (Retired)	Yes	Yes
Ultrasound Multi-frame Image Storage	Yes	Yes
Ultrasound Multi-frame Image Storage (Retired)	Yes	Yes
Video Endoscopic Image Storage	Yes	Yes
Video Microscopic Image Storage	Yes	Yes
Video Photographic Image Storage	Yes	Yes
VL Endoscopic Image Storage	Yes	Yes
VL Microscopic Image Storage	Yes	Yes
VL Photographic Image Storage	Yes	Yes
VL Slide-Coordinates Microscopic Image Storage	Yes	Yes
X-Ray 3D Angiographic Image Storage	Yes	Yes
X-Ray 3D Craniofacial Image Storage	Yes	Yes
X-Ray Angiographic Bi-Plane Image Storage (Retired)	Yes	Yes
X-Ray Angiographic Image Storage	Yes	Yes
X-Ray Radiofluoroscopic Image Storage	Yes	Yes
Query/Retrieve		
Study Root Information Model FIND	No	Yes
Study Root Information Model MOVE	No	Yes
Patient Root Information Model FIND	No	Yes
Patient Root Information Model MOVE	No	Yes
Waveforms, Notes, Reports, Measurements Transfer		
12-lead ECG Waveform Storage	Yes	Yes
Ambulatory ECG Waveform Storage	Yes	Yes
Basic Text SR	Yes	Yes
Basic Voice Audio Waveform Storage	Yes	Yes
Blending Softcopy Presentation State Storage SOP Class	Yes	Yes
Cardiac Electrophysiology Waveform Storage	Yes	Yes
Chest CAD SR	Yes	Yes
Color Softcopy Presentation State Storage SOP Class	Yes	Yes
Comprehensive SR	Yes	Yes
Deformable Spatial Registration Storage	Yes	Yes
Encapsulated CDA Storage	Yes	Yes
Encapsulated PDF Storage	Yes	Yes
Enhanced SR	Yes	Yes
General ECG Waveform Storage	Yes	Yes
Grayscale Softcopy Presentation State Storage SOP Class	Yes	Yes

Hemodynamic Waveform Storage	Yes	Yes
Key Object Selection Document	Yes	Yes
Mammography CAD SR	Yes	Yes
MR Spectroscopy Storage	Yes	Yes
Procedure Log Storage	Yes	Yes
Pseudo-Color Softcopy Presentation State Storage SOP Class	Yes	Yes
Raw Data Storage	Yes	Yes
Real World Value Mapping Storage	Yes	Yes
RT Beams Treatment Record Storage	Yes	Yes
RT Brachy Treatment Record Storage	Yes	Yes
RT Dose Storage	Yes	Yes
RT Ion Beams Treatment Record Storage	Yes	Yes
RT Ion Plan Storage	Yes	Yes
RT Plan Storage	Yes	Yes
RT Structure Set Storage	Yes	Yes
RT Treatment Summary Record Storage	Yes	Yes
Segmentation Storage	Yes	Yes
Spatial Fiducials Storage	Yes	Yes
Spatial Registration Storage	Yes	Yes
Standalone Curve Storage (Retired)	Yes	Yes
Standalone Modality LUT Storage (Retired)	Yes	Yes
Standalone Overlay Storage (Retired)	Yes	Yes
Standalone PET Curve Storage (Retired)	Yes	Yes
Standalone VOI LUT Storage (Retired)	Yes	Yes
Stereometric Relationship Storage	Yes	Yes
Stored Print Storage SOP Class (Retired)	Yes	Yes
X-Ray Radiation Dose SR	Yes	Yes

2. TABLE OF CONTENTS

Table of Contents

1. .. COMFORMANCE STATEMENT OVERVIEW	2
2. .. TABLE OF CONTENTS	5
3. .. INTRODUCTION	7
3.1 REVISION HISTORY	7
3.2 AUDIENCE	7
3.3 REMARKS	7
3.4 TERMS AND DEFINITIONS	7
3.5 BASICS OF DICOM COMMUNICATION	9
3.6 ABBREVIATIONS	9
3.7 REFERENCES	11
4. .. NETWORKING	12
4.1 IMPLEMENTATION MODEL	12
4.1.1Application Data Flow	12
4.1.2Functional Definitions of AE's	13
4.1.3Sequencing of Real-World Activities	13
4.2 AE SPECIFICATIONS	13
4.2.1ECHO-SCP	13
4.2.2STORAGE-SCP	15
4.2.3FIND-SCP	20
4.2.4MOVE-SCP	25
4.2.5STORAGE-SCU	28
4.3 NETWORK INTERFACES	33
4.3.1Physical Network Interface	33
4.3.2Additional Protocols	33
4.3.3IPv4 and IPv6 Support	34
4.4 CONFIGURATION	34
4.4.1AE Title/Presentation Address Mapping	34
4.4.2Parameters	34
5. .. MEDIA INTERCHANGE	35
6. .. SUPPORT OF CHARACTER SETS	36
6.1 OVERVIEW	36
6.2 CHARACTER SETS	36
6.3 CHARACTER SET CONFIGURATION	37
7. .. SECURITY	38
7.1 SECURITY PROFILES	38
7.2 ASSOCIATION LEVEL SECURITY	38
7.3 APPLICATION LEVEL SECURITY	38
8. .. ANNEXES	39
8.1 IOD CONTENTS	39
8.1.1Coerced/Modified fields	39
8.2 DATA DICTIONARY OF PRIVATE ATTRIBUTES	39
8.3 CODED TERMINOLOGY AND TEMPLATES	39
8.4 GRAYSCALE IMAGE CONSISTENCY	39
8.5 STANDARD EXTENDED/SPECIALIZED/PRIVATE SOP CLASSES	39
8.6 PRIVATE TRANSFER SYNTAXES	39

3. INTRODUCTION

3.1 REVISION HISTORY

Document Version	Date of Issue	Author	Description
1.0	December 11, 2008	Steve Wranovsky	Version for Final Review

3.2 AUDIENCE

This document is written for the people that need to understand how ClearCanvas ImageServer will integrate into their healthcare facility. This includes both those responsible for overall imaging network policy and architecture, as well as integrators who need to have a detailed understanding of the DICOM features of the product. This document contains some basic DICOM definitions so that any reader may understand how this product implements DICOM features. However, integrators are expected to fully understand all the DICOM terminology, how the tables in this document relate to the product's functionality, and how that functionality integrates with other devices that support compatible DICOM features.

3.3 REMARKS

The scope of this DICOM Conformance Statement is to facilitate integration between the ClearCanvas ImageServer and other DICOM products. The Conformance Statement should be read and understood in conjunction with the DICOM Standard. DICOM by itself does not guarantee interoperability. The Conformance Statement does, however, facilitate a first-level comparison for interoperability between different applications supporting compatible DICOM functionality.

This Conformance Statement is not supposed to replace validation with other DICOM equipment to ensure proper exchange of intended information. In fact, the user should be aware of the following important issues:

- The comparison of different Conformance Statements is just the first step towards assessing interconnectivity and interoperability between the product and other DICOM conformant equipment.
- Test procedures should be defined and executed to validate the required level of interoperability with specific compatible DICOM equipment, as established by the healthcare facility.

3.4 TERMS AND DEFINITIONS

Informal definitions are provided for the following terms used in this Conformance Statement. The DICOM Standard is the authoritative source for formal definitions of these terms.

Abstract Syntax – the information agreed to be exchanged between applications, generally equivalent to a Service/Object Pair (SOP) Class. Examples : Verification SOP Class, Modality Worklist Information Model Find SOP Class, Computed Radiography Image Storage SOP Class.

Application Entity (AE) – an end point of a DICOM information exchange, including the DICOM network or media interface software; i.e., the software that sends or receives DICOM information objects or messages. A single device may have multiple Application Entities.

Application Entity Title – the externally known name of an *Application Entity*, used to identify a DICOM application to other DICOM applications on the network.

Application Context – the specification of the type of communication used between *Application Entities*. Example: DICOM network protocol.

Association – a network communication channel set up between *Application Entities*.

Attribute – a unit of information in an object definition; a data element identified by a *tag*. The information may be a complex data structure (Sequence), itself composed of lower level data elements. Examples: Patient ID (0010,0020), Accession Number (0008,0050), Photometric Interpretation (0028,0004), Procedure Code Sequence (0008,1032).

Information Object Definition (IOD) – the specified set of *Attributes* that comprise a type of data object; does not represent a specific instance of the data object, but rather a class of similar data objects that have the same properties. The *Attributes* may be specified as Mandatory (Type 1), Required but possibly unknown (Type 2), or Optional (Type 3), and there may be conditions associated with the use of an Attribute (Types 1C and 2C). Examples: MR Image IOD, CT Image IOD, Print Job IOD.

Joint Photographic Experts Group (JPEG) – a set of standardized image compression techniques, available for use by DICOM applications.

Media Application Profile – the specification of DICOM information objects and encoding exchanged on removable media (e.g., CDs).

Module – a set of *Attributes* within an *Information Object Definition* that are logically related to each other. Example: Patient Module includes Patient Name, Patient ID, Patient Birth Date, and Patient Sex.

Negotiation – first phase of *Association* establishment that allows *Application Entities* to agree on the types of data to be exchanged and how that data will be encoded.

Presentation Context – the set of DICOM network services used over an *Association*, as negotiated between *Application Entities*; includes *Abstract Syntaxes* and *Transfer Syntaxes*.

Protocol Data Unit (PDU) – a packet (piece) of a DICOM message sent across the network. Devices must specify the maximum size packet they can receive for DICOM messages.

Security Profile – a set of mechanisms, such as encryption, user authentication, or digital signatures, used by an *Application Entity* to ensure confidentiality, integrity, and/or availability of exchanged DICOM data.

Service Class Provider (SCP) – role of an *Application Entity* that provides a DICOM network service; typically, a server that performs operations requested by another *Application Entity* (*Service Class User*). Examples: Picture Archiving and Communication System (image storage SCP, and image query/retrieve SCP), Radiology Information System (modality worklist SCP).

Service Class User (SCU) – role of an *Application Entity* that uses a DICOM network service; typically, a client. Examples: imaging modality (image storage SCU, and modality worklist SCU), imaging workstation (image query/retrieve SCU).

Service/Object Pair (SOP) Class – the specification of the network or media transfer (service) of a particular type of data (object); the fundamental unit of DICOM interoperability specification. Examples: Ultrasound Image Storage Service, Basic Grayscale Print Management.

Service/Object Pair (SOP) Instance – an information object; a specific occurrence of information exchanged in a *SOP Class*. Examples: a specific x-ray image.

Tag – a 32-bit identifier for a data element, represented as a pair of four digit hexadecimal numbers, the “group” and the “element”. If the “group” number is odd, the tag is for a private (manufacturer-specific) data element. Examples: (0010,0020) [Patient ID], (07FE,0010) [Pixel Data], (0019,0210) [private data element]

Transfer Syntax – the encoding used for exchange of DICOM information objects and messages. Examples: *JPEG* compressed (images), little endian explicit value representation.

Unique Identifier (UID) – a globally unique “dotted decimal” string that identifies a specific object or a class of objects; an ISO-8824 Object Identifier. Examples: Study Instance UID, SOP Class UID, SOP Instance UID.

Value Representation (VR) – the format type of an individual DICOM data element, such as text, an integer, a person’s name, or a code. DICOM information objects can be transmitted with either explicit identification of the type of each data element (Explicit VR), or without explicit identification (Implicit VR); with Implicit VR, the receiving application must use a DICOM data dictionary to look up the format of each data element.

3.5 BASICS OF DICOM COMMUNICATION

This section describes terminology used in this Conformance Statement for the non-specialist. The key terms used in the Conformance Statement are highlighted in *italics* below. This section is not a substitute for training about DICOM, and it makes many simplifications about the meanings of DICOM terms.

Two *Application Entities* (devices) that want to communicate with each other over a network using DICOM protocol must first agree on several things during an initial network “handshake”. One of the two devices must initiate an *Association* (a connection to the other device), and ask if specific services, information, and encoding can be supported by the other device (*Negotiation*).

DICOM specifies a number of network services and types of information objects, each of which is called an *Abstract Syntax* for the Negotiation. DICOM also specifies a variety of methods for encoding data, denoted *Transfer Syntaxes*. The Negotiation allows the initiating Application Entity to propose combinations of Abstract Syntax and Transfer Syntax to be used on the Association; these combinations are called *Presentation Contexts*. The receiving Application Entity accepts the Presentation Contexts it supports.

For each Presentation Context, the Association Negotiation also allows the devices to agree on *Roles* – which one is the *Service Class User* (SCU - client) and which is the *Service Class Provider* (SCP - server). Normally the device initiating the connection is the SCU, i.e., the client system calls the server, but not always.

The Association Negotiation finally enables exchange of maximum network packet (*PDU*) size, security information, and network service options (called *Extended Negotiation* information).

The Application Entities, having negotiated the Association parameters, may now commence exchanging data. Common data exchanges include queries for worklists and lists of stored images, transfer of image objects and analyses (structured reports), and sending images to film printers. Each exchangeable unit of data is formatted by the sender in accordance with the appropriate *Information Object Definition*, and sent using the negotiated Transfer Syntax. There is a Default Transfer Syntax that all systems must accept, but it may not be the most efficient for some use cases. Each transfer is explicitly acknowledged by the receiver with a *Response Status* indicating success, failure, or that query or retrieve operations are still in process.

Two Application Entities may also communicate with each other by exchanging media (such as a CD-R). Since there is no Association Negotiation possible, they both use a *Media Application Profile* that specifies “pre-negotiated” exchange media format, Abstract Syntax, and Transfer Syntax.

3.6 ABBREVIATIONS

Abbreviations should be listed here. These may be taken from the following list, deleting terms that are not used within the Conformance Statement, and adding any additional terms that are used:

AE Application Entity

AET	Application Entity Title
CAD	Computer Aided Detection
CDA	Clinical Document Architecture
CD-R	Compact Disk Recordable
CSE	Customer Service Engineer
CR	Computed Radiography
CT	Computed Tomography
DHCP	Dynamic Host Configuration Protocol
DICOM	Digital Imaging and Communications in Medicine
DNS	Domain Name System
DX	Digital X-ray
GSPTS	Grayscale Softcopy Presentation State
HIS	Hospital Information System
HL7	Health Level 7 Standard
IHE	Integrating the Healthcare Enterprise
IOD	Information Object Definition
IPv4	Internet Protocol version 4
IPv6	Internet Protocol version 6
ISO	International Organization for Standards
IO	Intra-oral X-ray
JPEG	Joint Photographic Experts Group
LUT	Look-up Table
MPEG	Moving Picture Experts Group
MG	Mammography (X-ray)
MPPS	Modality Performed Procedure Step
MR	Magnetic Resonance Imaging
MSPS	Modality Scheduled Procedure Step
MTU	Maximum Transmission Unit (IP)
MWL	Modality Worklist
NM	Nuclear Medicine
NTP	Network Time Protocol
O	Optional (Key Attribute)
OP	Ophthalmic Photography
OSI	Open Systems Interconnection
PACS	Picture Archiving and Communication System
PET	Positron Emission Tomography
PDU	Protocol Data Unit
R	Required (Key Attribute)
RF	Radiofluoroscopy
RIS	Radiology Information System.
RT	Radiotherapy
SC	Secondary Capture
SCP	Service Class Provider
SCU	Service Class User

SOP	Service-Object Pair
SPS	Scheduled Procedure Step
SR	Structured Reporting
TCP/IP	Transmission Control Protocol/Internet Protocol
U	Unique (Key Attribute)
UL	Upper Layer
US	Ultrasound
VL	Visible Light
VR	Value Representation
XA	X-ray Angiography


3.7 REFERENCES

- NEMA PS3 Digital Imaging and Communications in Medicine (DICOM) Standard, available free at <http://medical.nema.org/>

4. NETWORKING

4.1 IMPLEMENTATION MODEL

4.1.1 Application Data Flow



**Figure 4.1-1
IMPLEMENTATION MODEL**

The application is a .NET application that provides a Web based user interface, internal database and network listener that spawns additional threads as necessary to handle incoming connections.

Conceptually, the network services may be modeled as the following separate AEs, though in fact all the AEs share a single (configurable) AE Title:

- ECHO-SCP, which responds to verification requests
- STORAGE-SCP, which receives incoming composite instances
- FIND-SCP, which receives incoming queries for lists of studies

- MOVE-SCP, which responds to requests for studies
- STORAGE-SCU, which sends outbound composite instances

4.1.2 Functional Definitions of AE’s

4.1.2.1 ECHO-SCP

ECHO-SCP waits in the background for connections, will accept associations with Presentation Contexts for the SOP Class of the Verification Service Class, and will respond successfully to echo requests.

4.1.2.2 STORAGE-SCP

STORAGE-SCP waits in the background for connections, will accept associations with Presentation Contexts for SOP Classes of the Storage Service Class, and will store the received instances to the local database where they may subsequently be listed and viewed through the user interface.

4.1.2.3 FIND-SCP

FIND- SCP waits in the background for connections, will accept associations with Presentation Contexts for the SOP Class of the Study Root Query/Retrieve Information Model – FIND Service Class or the Patient Root Query/Retrieve Information Model – FIND Service Class, and will respond successfully to query requests.

4.1.2.4 MOVE-SCP

MOVE-SCP waits in the background for connections, will accept associations with Presentation Contexts for the SOP Class of the Study Root Query/Retrieve Information Model – MOVE Service Class or the Patient Root Query/Retrieve Information Model – MOVE Service Class, and will respond successfully to retrieve requests by initiating storage of instances to the remote Application Entity.

4.1.2.5 STORAGE-SCU

STORAGE-SCU is activated through the user interface when a user selects studies from the local database and requests that they be sent to a remote AE (selected from a pre-configured list).

4.1.3 Sequencing of Real-World Activities

All SCP activities are performed asynchronously in the background and are not dependent on any sequencing.

All SCU activities are initiated through the user interface, with the exception of STORAGE-SCU which is also initiated by MOVE-SCP in order to store the requested instances.

4.2 AE SPECIFICATIONS

4.2.1 ECHO-SCP

4.2.1.1 SOP Classes

ECHO-SCP provides Standard Conformance to the following SOP Class(es):

**Table 4.2-1
SOP CLASSES SUPPORTED BY ECHO-SCP**

SOP Class Name	SOP Class UID
Verification SOP Class	1.2.840.10008.1.1

4.2.1.2 Association Policies

4.2.1.2.1 General

ECHO-SCP accepts but never initiates associations.

**Table 4.2-2
MAXIMUM PDU SIZE RECEIVED AS A SCP FOR ECHO-SCP**

Maximum PDU size received	114kB (approximate)
---------------------------	---------------------

4.2.1.2.2 Number of Associations

**Table 4.2-3
NUMBER OF ASSOCIATIONS AS A SCP FOR ECHO-SCP**

Maximum number of simultaneous associations	Unlimited
---	-----------

4.2.1.2.3 Asynchronous Nature

ECHO-SCP will only allow a single outstanding operation on an Association. Therefore, ECHO-SCP will not perform asynchronous operations window negotiation.

4.2.1.2.4 Implementation Identifying Information

**Table 4.2-4
DICOM IMPLEMENTATION CLASS AND VERSION FOR ECHO-SCP**

Implementation Class UID	1.3.6.1.4.1.25403.1.1.1
Implementation Version Name	Dicom 0.1

4.2.1.3 Association Initiation Policy

ECHO-SCP does not initiate associations.

4.2.1.4 Association Acceptance Policy

When ECHO-SCP accepts an association, it will respond to echo requests. If the Called AE Title does not match the pre-configured AE Title shared by all the SCPs of the application, the association will be rejected.

4.2.1.4.1 Activity – Receive Echo Request

4.2.1.4.1.1 Description and Sequencing of Activities

As requests are received, they are responded to immediately.

4.2.1.4.1.2 Accepted Presentation Contexts

**Table 4.2-5
ACCEPTABLE PRESENTATION CONTEXTS FOR ECHO-SCP AND RECEIVE ECHO REQUEST**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID		
Verification	1.2.840.10008.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None

4.2.1.4.1.2.1 Extended Negotiation

No extended negotiation is performed.

4.2.1.4.1.3 SOP Specific Conformance

4.2.1.4.1.3.1 SOP Specific Conformance to Verification SOP Class

ECHO-SCP provides standard conformance to the Verification Service Class.

4.2.1.4.1.3.2 Presentation Context Acceptance Criterion

ECHO-SCP will only accept a Presentation Context compatible with the one listed in Table 4.2-5.

4.2.1.4.1.3.3 Transfer Syntax Selection Policies

If proposed, ECHO-SCP prefers the Explicit VR Little Endian Transfer Syntax.

ECHO-SCP will accept duplicate Presentation Contexts; that is, if it is offered multiple Presentation Contexts, each of which offers acceptable Transfer Syntaxes, it will accept all Presentation Contexts, applying the same method for selecting a Transfer Syntax for each.

4.2.2 STORAGE-SCP

4.2.2.1 SOP Classes

STORAGE-SCP provides Standard Conformance to the following SOP Class(es):

**Table 4.2-6
IMAGE SOP CLASSES SUPPORTED BY STORAGE-SCP**

SOP Class Name	SOP Class UID
Computed Radiography Image Storage	1.2.840.10008.5.1.4.1.1.1
CT Image Storage	1.2.840.10008.5.1.4.1.1.2
Digital Intra-oral X-Ray Image Storage – For Presentation	1.2.840.10008.5.1.4.1.1.1.3
Digital Intra-oral X-Ray Image Storage – For Processing	1.2.840.10008.5.1.4.1.1.1.3.1
Digital Mammography X-Ray Image Storage – For Presentation	1.2.840.10008.5.1.4.1.1.1.2
Digital Mammography X-Ray Image Storage – For Processing	1.2.840.10008.5.1.4.1.1.1.2.1
Digital X-Ray Image Storage – For Presentation	1.2.840.10008.5.1.4.1.1.1.1
Digital X-Ray Image Storage – For Processing	1.2.840.10008.5.1.4.1.1.1.1.1
Enhanced CT Image Storage	1.2.840.10008.5.1.4.1.1.2.1
Enhanced MR Image Storage	1.2.840.10008.5.1.4.1.1.4.1
Enhanced XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1.1
Enhanced XRF Image Storage	1.2.840.10008.5.1.4.1.1.12.2.1
Hardcopy Grayscale Image Storage SOP Class (Retired)	1.2.840.10008.5.1.1.29
Hardcopy Color Image Storage SOP Class (Retired)	1.2.840.10008.5.1.1.30
MR Image Storage	1.2.840.10008.5.1.4.1.1.4
Multi-frame 1 Color Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.4
Multi-frame Grayscale Byte Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.2
Multi-frame Grayscale Word Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.3
Multi-frame Single Bit Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.1
Nuclear Medicine Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.5
Nuclear Medicine Image Storage	1.2.840.10008.5.1.4.1.1.20

Ophthalmic Photography 16 Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.2
Ophthalmic Photography 8 Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.1
Ophthalmic Tomography Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.4
Positron Emission Tomography Image Storage	1.2.840.10008.5.1.4.1.1.128
RT Image Storage	1.2.840.10008.5.1.4.1.1.481.1
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1
Ultrasound Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.6
Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1
Ultrasound Multi-frame Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.3
Video Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1.1
Video Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2.1
Video Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4.1
VL Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1
VL Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2
VL Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4
VL Slide-Coordinates Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.3
X-Ray 3D Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.13.1.1
X-Ray 3D Craniofacial Image Storage	1.2.840.10008.5.1.4.1.1.13.1.2
X-Ray Angiographic Bi-Plane Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.12.3
X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1
X-Ray Radiofluoroscopic Image Storage	1.2.840.10008.5.1.4.1.1.12.2

**Table 4.2-7
NON-IMAGE SOP CLASSES SUPPORTED BY STORAGE-SCP**

SOP Class Name	SOP Class UID
12-lead ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1
Ambulatory ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.3
Basic Text SR	1.2.840.10008.5.1.4.1.1.88.11
Basic Voice Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.1
Blending Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.4
Cardiac Electrophysiology Waveform Storage	1.2.840.10008.5.1.4.1.1.9.3.1
Chest CAD SR	1.2.840.10008.5.1.4.1.1.88.65
Color Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.2
Comprehensive SR	1.2.840.10008.5.1.4.1.1.88.33
Deformable Spatial Registration Storage	1.2.840.10008.5.1.4.1.1.66.3
Encapsulated CDA Storage	1.2.840.10008.5.1.4.1.1.104.2
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1
Enhanced SR	1.2.840.10008.5.1.4.1.1.88.22

General ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.2
Grayscale Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.1
Hemodynamic Waveform Storage	1.2.840.10008.5.1.4.1.1.9.2.1
Key Object Selection Document	1.2.840.10008.5.1.4.1.1.88.59
Mammography CAD SR	1.2.840.10008.5.1.4.1.1.88.50
MR Spectroscopy Storage	1.2.840.10008.5.1.4.1.1.4.2
Procedure Log Storage	1.2.840.10008.5.1.4.1.1.88.40
Pseudo-Color Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.3
Raw Data Storage	1.2.840.10008.5.1.4.1.1.66
Real World Value Mapping Storage	1.2.840.10008.5.1.4.1.1.67
RT Beams Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.4
RT Brachy Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.6
RT Dose Storage	1.2.840.10008.5.1.4.1.1.481.2
RT Ion Beams Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.9
RT Ion Plan Storage	1.2.840.10008.5.1.4.1.1.481.8
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5
RT Structure Set Storage	1.2.840.10008.5.1.4.1.1.481.3
RT Treatment Summary Record Storage	1.2.840.10008.5.1.4.1.1.481.7
Segmentation Storage	1.2.840.10008.5.1.4.1.1.66.4
Spatial Fiducials Storage	1.2.840.10008.5.1.4.1.1.66.2
Spatial Registration Storage	1.2.840.10008.5.1.4.1.1.66.1
Standalone Curve Storage (Retired)	1.2.840.10008.5.1.4.1.1.9
Standalone Modality LUT Storage (Retired)	1.2.840.10008.5.1.4.1.1.10
Standalone Overlay Storage (Retired)	1.2.840.10008.5.1.4.1.1.8
Standalone PET Curve Storage (Retired)	1.2.840.10008.5.1.4.1.1.129
Standalone VOI LUT Storage (Retired)	1.2.840.10008.5.1.4.1.1.11
Stereometric Relationship Storage	1.2.840.10008.5.1.4.1.1.77.1.5.3
Stored Print Storage SOP Class (Retired)	1.2.840.10008.5.1.1.27
X-Ray Radiation Dose SR	1.2.840.10008.5.1.4.1.1.88.67

4.2.2.2 Association Policies

4.2.2.2.1 General

STORAGE-SCP accepts but never initiates associations.

**Table 4.2-8
MAXIMUM PDU SIZE RECEIVED AS A SCP FOR STORAGE-SCP**

Maximum PDU size received	114kB (approx.)
---------------------------	-----------------

4.2.2.2 Number of Associations

**Table 4.2-9
NUMBER OF ASSOCIATIONS AS A SCP FOR STORAGE-SCP**

Maximum number of simultaneous associations	Unlimited
---	-----------

4.2.2.2.3 Asynchronous Nature

STORAGE-SCP will only allow a single outstanding operation on an Association. Therefore, STORAGE-SCP will not perform asynchronous operations window negotiation.

4.2.2.2.4 Implementation Identifying Information

**Table 4.2-10
DICOM IMPLEMENTATION CLASS AND VERSION FOR STORAGE-SCP**

Implementation Class UID	1.3.6.1.4.1.25403.1.1.1
Implementation Version Name	Dicom 0.1

4.2.2.3 Association Initiation Policy

STORAGE-SCP does not initiate associations.

4.2.2.4 Association Acceptance Policy

When STORAGE-SCP accepts an association, it will respond to storage requests. If the Called AE Title does not match the pre-configured AE Title shared by all the SCPs of the application, the association will be rejected.

4.2.2.4.1 Activity – Receive Storage Request

4.2.2.4.2 Description and Sequencing of Activities

As instances are received, they are copied to the local file system and a record inserted into the local database. The STORAGE-SCP is configurable as to how it handles duplicate SOP Instances. It can accept, reject, ignore, or compare the SOP instances to the duplicates.

4.2.2.4.2.1 Accepted Presentation Contexts

Table 4.2-11 contains the transfer syntaxes supported for Image SOP Classes and Table 4.2-12 contains the transfer syntaxes supported for Non-Image SOP Classes.

**Table 4.2-11
ACCEPTABLE PRESENTATION CONTEXTS FOR
STORAGE-SCP AND RECEIVE STORAGE REQUEST FOR IMAGE SOP CLASSES**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID		
See Table 4.2-6	See Table 4.2-6	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None
		RLE Lossless	1.2.840.10008.1.2.5	SCP	None
		JPEG Lossless, Non-Hierarchical (Process 14)	1.2.840.10008.1.2.4.57	SCP	None

	JPEG Extended (Process 2 & 4)	1.2.840.10008.1.2.4.51	SCP	None
	JPEG Baseline (Process 1):	1.2.840.10008.1.2.4.50	SCP	None
	JPEG 2000 Image Compression (Lossless Only)	1.2.840.10008.1.2.4.90	SCP	None
	JPEG 2000 Image Compression	1.2.840.10008.1.2.4.91	SCP	None

**Table 4.2-12
ACCEPTABLE PRESENTATION CONTEXTS FOR
STORAGE-SCP AND RECEIVE STORAGE REQUEST FOR NON-IMAGE SOP CLASSES**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID		
See Table 4.2-7	See Table 4.2-7	Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None
		Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None

4.2.2.4.2.1.1 Extended Negotiation

No extended negotiation is performed, though STORAGE-SCP:

- is a Level 2 Storage SCP (Full – does not discard any data elements)
- does not support digital signatures
- In some cases, it may coerce the Patient Name, Patient ID, Issuer of Patient ID, Patient Birth Date, Patient Sex, and Accession Number fields. See section 8.1.1 for details.

4.2.2.4.2.2 SOP Specific Conformance

4.2.2.4.2.2.1 SOP Specific Conformance to Storage SOP Classes

STORAGE-SCP provides standard conformance to the Storage Service Class.

4.2.2.4.2.2.2 Presentation Context Acceptance Criterion

STORAGE-SCP will always accept any Presentation Context for the supported SOP Classes with the supported Transfer Syntaxes. More than one proposed Presentation Context will be accepted for the same Abstract Syntax if the Transfer Syntax is supported, whether or not it is the same as another Presentation Context.

4.2.2.4.2.2.3 Transfer Syntax Selection Policies

The STORAGE-SCP will always accept the first acceptable presentation context proposed by the SCU.

STORAGE-SCP will accept duplicate Presentation Contexts; that is, if it is offered multiple Presentation Contexts, each of which offers acceptable Transfer Syntaxes, it will accept all Presentation Contexts, applying the same method for selecting a Transfer Syntax for each.

4.2.2.4.2.4 Response Status

STORAGE-SCP will behave as described in the Table below when generating the C-STORE response command message.

**Table 4.2-13
RESPONSE STATUS FOR STORAGE-SCP AND RECEIVE STORAGE REQUEST**

Service Status	Further Meaning	Status Codes	Reason
Failure	Refused: Out of Resources	A700	STORAGE-SCP does not have an online filesystem with sufficient space to store the incoming instance or the Study is currently being processed and cannot accept new instances at this time, or the study has been lossy compressed and cannot accept new instances.
	Duplicate SOP Instance	0111	The SOP Instance has been received already and the STORAGE-SCP is configured to reject duplicates.
	Attribute Value Out of Range	0116	SOP Instance UID, Study Instance UID, or Series Instance UID are longer than allowed by DICOM (greater than 64 characters)
	Processing Failure	0110	A failure was encountered when processing the SOP Instance.
Success		0000	

4.2.3 FIND-SCP

4.2.3.1 SOP Classes

FIND-SCP provides Standard Conformance to the following SOP Class(es):

**Table 4.2-14
SOP CLASSES SUPPORTED BY FIND-SCP**

SOP Class Name	SOP Class UID
Study Root Query/Retrieve Information Model – FIND	1.2.840.10008.5.1.4.1.2.2.1
Patient Root Query/Retrieve Information Model – FIND	1.2.840.10008.5.1.4.1.2.1.1

4.2.4.2 Association Policies

4.2.3.1.1 General

FIND-SCP accepts but never initiates associations.

**Table 4.2-15
Maximum PDU size received as a SCP for FIND-SCP**

Maximum PDU size received	114kB (approx.)
---------------------------	-----------------

4.2.3.1.2 Number of Associations

**Table 4.2-16
Number of Associations as a SCP for FIND-SCP**

Maximum number of simultaneous associations	Unlimited
---	-----------

4.2.3.1.3 Asynchronous Nature

FIND-SCP will only allow a single outstanding operation on an Association. Therefore, FIND-SCP will not perform asynchronous operations window negotiation.

4.2.3.1.4 Implementation Identifying Information

**Table 4.2-17
DICOM Implementation Class and Version for FIND-SCP**

Implementation Class UID	1.3.6.1.4.1.25403.1.1.1
Implementation Version Name	Dicom 0.1

4.2.3.2 Association Initiation Policy

FIND-SCP does not initiate associations.

4.2.3.3 Association Acceptance Policy

When FIND-SCP accepts an association, it will respond to query requests. If the Called AE Title does not match the pre-configured AE Title shared by all the SCPs of the application, the association will be rejected.

4.2.3.3.1 Activity – Receive Query Request

4.2.3.3.1.1 Description and Sequencing of Activities

When a query is received, the local database is queried for the result set.

4.2.3.3.1.2 Accepted Presentation Contexts

**Table 4.2-18
Acceptable Presentation Contexts for FIND-SCP and Incoming Query from Remote AE**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID		
See Table 4.2-14	See Table 4.2-14	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None

4.2.3.3.1.2.1 Extended Negotiation

No extended negotiation is performed.

In particular, relational queries are not supported.

4.2.3.3.1.3 SOP Specific Conformance

4.2.3.3.1.3.1 SOP Specific Conformance to C-FIND SOP Classes

FIND-SCP provides standard conformance to the supported C-FIND SOP Classes.

Only those attributes that are requested are returned in a C-FIND response. Some optional requested attributes will be returned for Patient Root queries as per Table 4.2-19.

**Table 4.2-19
PATIENT ROOT RESPONSE IDENTIFIER FOR FIND-SCP**

Name	Tag	Types of Matching
PATIENT Level		
Patient's ID	(0010,0020)	S,*,U
Patient's Name	(0010,0010)	S,*,U
Issuer of Patient ID	(0010,0021)	S,*,U
STUDY Level		
Study Instance UID	(0020,000D)	UNIQUE
Study ID	(0020,0010)	S,*,U
Study Description	(0008,1030)	S,*,U
Modalities In Study	(0008,0061)	S,*,U
Study Date	(0008,0020)	S,U,R
Study Time	(0008,0030)	R,U
Accession Number	(0008,0050)	S,*,U
Referring Physician's Name	(0008,0090)	S,*,U
Number Of Study Related Series	(0020, 1206)	S,*,U
Number Of Study Related Instances	(0020, 1208)	S,*,U
SERIES Level		
Series Instance UID	(0020, 000E)	UNIQUE
Modality	(0008, 0060)	S,*,U
Series Description	(0008, 103E)	S,*,U
Series Number	(0020, 0011)	S,*,U
Performed Procedure Step Start Date	(0040,0244)	R, U
Performed Procedure Step Start Time	(0040,0245)	R, U
Request Attributes Sequence	(0040,0275)	NONE
>Scheduled Procedure Step ID	(0040,0009)	NONE
>Requested Procedure ID	(0040,1001)	NONE
Number Of Series Related Instances	(0020, 1209)	S,*,U
IMAGE Level		
Sop Instance UID	(0008, 0018)	UNIQUE
The FIND-SCP supports queries for any Image Level DICOM attribute with the exception of Pixel Data.		
Common to all query levels		

Specific Character Set	(0008,0005)	N/A
------------------------	-------------	-----

Only those attributes that are requested are returned in a C-FIND response. Some optional requested attributes will be returned for Study Root queries as per Table 4.2-20.

**Table 4.2-20
STUDY ROOT RESPONSE IDENTIFIER FOR FIND-SCP**

Name	Tag	Types of Matching
STUDY Level		
Study Instance UID	(0020,000D)	UNIQUE
Patient's ID	(0010,0020)	S,*,U
Patient's Name	(0010,0010)	S,*,U
Patient's Birth Date	(0010,0030)	R,U
Patient's Sex	(0010,0040)	S,*,U
Study ID	(0020,0010)	S,*,U
Study Description	(0008,1030)	S,*,U
Modalities In Study	(0008,0061)	S,*,U
Study Date	(0008,0020)	S,U,R
Study Time	(0008,0030)	R,U
Accession Number	(0008,0050)	S,*,U
Referring Physician's Name	(0008,0090)	S,*,U
Number Of Study Related Series	(0020, 1206)	S,*,U
Number Of Study Related Instances	(0020, 1208)	S,*,U
SERIES Level		
Series Instance UID	(0020, 000E)	UNIQUE
Modality	(0008, 0060)	S,*,U
Series Description	(0008, 103E)	S,*,U
Series Number	(0020, 0011)	S,*,U
Performed Procedure Step Start Date	(0040,0244)	R, U
Performed Procedure Step Start Time	(0040,0245)	R, U
Request Attributes Sequence	(0040,0275)	NONE
>Scheduled Procedure Step ID	(0040,0009)	NONE
>Requested Procedure ID	(0040,1001)	NONE
Number Of Series Related Instances	(0020, 1209)	S,*,U
IMAGE Level		
Sop Instance UID	(0008, 0018)	UNIQUE
The FIND-SCP supports queries for any Image Level DICOM attribute with the exception of Pixel Data.		
Common to all query levels		
Specific Character Set	(0008,0005)	N/A

Types of Matching:

An "S" indicates the identifier attribute uses Single Value Matching, an "R" indicates Range Matching, a "*" indicates wildcard matching, a 'U' indicates Universal Matching, and an 'L' indicates that UID lists are sent. "NONE" indicates that no matching is supported, but that values for this Element are requested to be returned (i.e. universal matching), and "UNIQUE" indicates that this is the Unique Key for that query level, in which case Universal Matching or Single Value Matching is used depending on the query level.

4.2.3.3.1.3.2 Presentation Context Acceptance Criterion

FIND-SCP will always accept any Presentation Context for the supported SOP Classes with the supported Transfer Syntaxes. More than one proposed Presentation Context will be accepted for the same Abstract Syntax if the Transfer Syntax is supported, whether or not it is the same as another Presentation Context.

4.2.3.3.1.3.3 Transfer Syntax Selection Policies

FIND-SCP will always accept the first acceptable proposed transfer syntax.

FIND-SCP will accept duplicate Presentation Contexts; that is, if it is offered multiple Presentation Contexts, each of which offers acceptable Transfer Syntaxes, it will accept all Presentation Contexts, applying the same method for selecting a Transfer Syntax for each.

4.2.3.3.1.3.4 Response Status

FIND-SCP will behave as described in Table 4.2-29 in response to the status returned in the C-FIND response command message(s).

**Table 4.2-21
RESPONSE STATUS FOR FIND-SCP AND RECEIVE QUERY REQUEST**

Service Status	Further Meaning	Status Codes	Behavior
Failure	Unable to process	C000	Sent if internal database query is unsuccessful
	Identifier does not match SOP Class	A900	Sent if an invalid query retrieve level is set.
Failure	Unable to process	C000	Sent if internal database query is unsuccessful
Cancel	Matching terminated due to Cancel request	FE00	Never sent. Cancel requests are ignored.
Success	Matching is complete	0000	Sent
Pending	Matches are continuing - Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys	FF00	Sent

4.2.4 MOVE-SCP

4.2.4.1 SOP Classes

MOVE-SCP provides Standard Conformance to the following SOP Class(es):

Table 4.2-22
SOP CLASSES SUPPORTED BY MOVE-SCP

SOP Class Name	SOP Class UID
Study Root Query/Retrieve Information Model – MOVE	1.2.840.10008.5.1.4.1.2.2.2
Patient Root Query/Retrieve INFORMATION Model -- MOVE	1.2.840.10008.5.1.4.1.2.1.2

4.2.4.2 Association Policies

4.2.4.2.1 General

MOVE-SCP accepts but never initiates associations.

Table 4.2-23
MAXIMUM PDU SIZE RECEIVED AS A SCP FOR MOVE-SCP

Maximum PDU size received	114kB (approx.)
---------------------------	-----------------

4.2.4.2.2 Number of Associations

Table 4.2-24
NUMBER OF ASSOCIATIONS AS A SCP FOR MOVE-SCP

Maximum number of simultaneous associations	Unlimited
---	-----------

4.2.4.2.3 Asynchronous Nature

MOVE-SCP will only allow a single outstanding operation on an Association. Therefore, MOVE-SCP will not perform asynchronous operations window negotiation.

4.2.4.2.4 Implementation Identifying Information

Table 4.2-25
DICOM IMPLEMENTATION CLASS AND VERSION FOR MOVE-SCP

Implementation Class UID	1.3.6.1.4.1.25403.1.1.1
Implementation Version Name	Dicom 0.1

4.2.4.3 Association Initiation Policy

MOVE-SCP initiates an association with the AE specified as the Move Destination in the MOVE request, in order to store the requested instances. The remote AE must be in the application's pre-configured AE list.

4.2.4.4 Association Acceptance Policy

When MOVE-SCP accepts an association, it will respond to move requests. If configured, the association will be rejected if the Calling AE Title is not pre-configured in the database.

4.2.4.4.1 Activity – Receive Move Request

4.2.4.4.1.1 Description and Sequencing of Activities

As requests are received, a STORAGE-SCU operation is initiated to send the requested instances to the specified remote AE.

4.2.4.4.1.2 Proposed Presentation Contexts

**Table 4.2-26
ACCEPTABLE PRESENTATION CONTEXTS FOR MOVE-SCP AND RECEIVE RETRIEVE REQUEST**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID		
See Table 4.2-22	See Table 4.2-22	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None

4.2.4.4.1.2.1 Extended Negotiation

No extended negotiation is performed through MOVE-SCP.

In particular, relational retrievals are not supported.

4.2.4.4.1.3 SOP Specific Conformance

4.2.4.4.1.3.1 SOP Specific Conformance to C-MOVE SOP Classes

MOVE-SCP provides standard conformance to the supported C-MOVE SOP Classes.

The Patient Root and Study Root information models are supported.

The move is performed to the destination AE Title specified in the original request. . If the destination AE does not exist in the application’s pre-configured list, the store operations are not performed.

**Table 4.2-27
PATIENT ROOT REQUEST IDENTIFIER FOR MOVE-SCP**

Name	Tag	Unique, Matching or Return Key
PATIENT level		
Patient ID	(0010,0020)	U
STUDY level		
Study Instance UID	(0020,000D)	U
SERIES level		
Study Instance UID	(0020,000D)	U
Series Instance UID	(0020,000E)	U
IMAGE level		
Study Instance UID	(0020,000D)	U
Series Instance UID	(0020,000E)	U

Sop Instance UID	(0008,0018)	U
------------------	-------------	---

**Table 4.2-28
STUDY ROOT REQUEST IDENTIFIER FOR MOVE-SCP**

Name	Tag	Unique, Matching or Return Key
STUDY level		
Study Instance UID	(0020,000D)	U
SERIES level		
Study Instance UID	(0020,000D)	U
Series Instance UID	(0020,000E)	U
IMAGE level		
Study Instance UID	(0020,000D)	U
Series Instance UID	(0020,000E)	U
Sop Instance UID	(0008,0018)	U

4.2.4.4.1.3.2 Presentation Context Acceptance Criterion

MOVE-SCP will only accept a Presentation Context compatible with the one listed in Table 4.2-26.

4.2.4.4.1.3.3 Transfer Syntax Selection Policies

MOVE -SCP will always select the first Transfer Syntax proposed by the client that is supported by the SCP, per Presentation Context.

MOVE -SCP will accept duplicate Presentation Contexts; that is, if it is offered multiple Presentation Contexts, each of which offers acceptable Transfer Syntaxes, it will accept all Presentation Contexts, applying the same method for selecting a Transfer Syntax for each.

4.2.4.4.1.3.4 Response Status

MOVE-SCP will behave as described in the Table below when generating the C-MOVE response command message.

**Table 4.2-29
RESPONSE STATUS FOR MOVE-SCP AND SEND TO REMOTE AE REQUEST**

Service Status	Further Meaning	Status Codes	Related Fields	Behavior
Failure	Refused: Out of Resources – Unable to perform sub-operations	A702	(0000,1020) (0000,1021) (0000,1022) (0000,1023)	The selected retrieve objects are not online, and cannot be retrieved at this time.
	Refused: Move Destination unknown	A801	(0000,0902)	Sent if the destination AE Title has not been preconfigured, or a connection failure occurs when connecting to the destination AE.

	Identifier does not match SOP Class	A900	(0000,0901) (0000,0902)	Sent when the Query/Retrieve Level is not present or is invalid
	Unable to process	C000	(0000,0901) (0000,0902)	Sent if the local database query fails or an internal exception occurred during processing
	Sub-operations terminated due to Cancel Indication	FE00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)	Sent if Move association is cancelled Sub-operations are also cancelled
	Sub-operations Complete - One or more Failures	B000	(0000,1020) (0000,1022) (0000,1023)	Sent
Success	Sub-operations Complete - No Failures	0000	(0000,1020) (0000,1021) (0000,1022) (0000,1023)	Sent
Pending	Sub-operations are continuing	FF00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)	Sent for every 5 C-STORE operations completed by STORAGE-SCU.

4.2.4.4.1.3.5 Sub-operation dependent behavior

Since the C-MOVE operation is dependent on completion of C-STORE sub-operations that are occurring on a separate association, the question of failure of operations on the other association(s) must be considered.

MOVE-SCP initiates a C-STORE sub-operation for each request. The responses from the MOVE-SCP are purely dependent on the success or failure of the C-STORE sub-operations, not on any explicit action by MOVE-SCP.

Whether or not the remote AE accepts the C-STORE sub-operations is beyond the control of MOVE-SCU.

If the association on which the C-MOVE was issued is aborted for any reason, the C-STORE sub-operations will continue.

If the C-MOVE operation is canceled by the remote AE, MOVE-SCP will also attempt to cancel the corresponding C-STORE sub-operation.

4.2.5 STORAGE-SCU

4.2.5.1 SOP Classes

STORAGE-SCU provide Standard Conformance to the following SOP Class(es):

**Table 4.2-30
IMAGE SOP CLASSES SUPPORTED BY STORAGE-SCU**

SOP Class Name	SOP Class UID
Computed Radiography Image Storage	1.2.840.10008.5.1.4.1.1.1

CT Image Storage	1.2.840.10008.5.1.4.1.1.2
Digital Intra-oral X-Ray Image Storage – For Presentation	1.2.840.10008.5.1.4.1.1.1.3
Digital Intra-oral X-Ray Image Storage – For Processing	1.2.840.10008.5.1.4.1.1.1.3.1
Digital Mammography X-Ray Image Storage – For Presentation	1.2.840.10008.5.1.4.1.1.1.2
Digital Mammography X-Ray Image Storage – For Processing	1.2.840.10008.5.1.4.1.1.1.2.1
Digital X-Ray Image Storage – For Presentation	1.2.840.10008.5.1.4.1.1.1.1
Digital X-Ray Image Storage – For Processing	1.2.840.10008.5.1.4.1.1.1.1.1
Enhanced CT Image Storage	1.2.840.10008.5.1.4.1.1.2.1
Enhanced MR Image Storage	1.2.840.10008.5.1.4.1.1.4.1
Enhanced XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1.1
Enhanced XRF Image Storage	1.2.840.10008.5.1.4.1.1.12.2.1
Hardcopy Grayscale Image Storage SOP Class (Retired)	1.2.840.10008.5.1.1.29
Hardcopy Color Image Storage SOP Class (Retired)	1.2.840.10008.5.1.1.30
MR Image Storage	1.2.840.10008.5.1.4.1.1.4
Multi-frame 1 Color Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.4
Multi-frame Grayscale Byte Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.2
Multi-frame Grayscale Word Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.3
Multi-frame Single Bit Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.1
Nuclear Medicine Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.5
Nuclear Medicine Image Storage	1.2.840.10008.5.1.4.1.1.20
Ophthalmic Photography 16 Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.2
Ophthalmic Photography 8 Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.1
Ophthalmic Tomography Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.4
Positron Emission Tomography Image Storage	1.2.840.10008.5.1.4.1.1.128
RT Image Storage	1.2.840.10008.5.1.4.1.1.481.1
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1
Ultrasound Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.6
Ultrasound Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1
Ultrasound Multi-frame Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.3
Video Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1.1
Video Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2.1
Video Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4.1
VL Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1
VL Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2
VL Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4
VL Slide-Coordinates Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.3
X-Ray 3D Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.13.1.1
X-Ray 3D Craniofacial Image Storage	1.2.840.10008.5.1.4.1.1.13.1.2

X-Ray Angiographic Bi-Plane Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.12.3
X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1
X-Ray Radiofluoroscopic Image Storage	1.2.840.10008.5.1.4.1.1.12.2

**Table 4.2-31
NON-IMAGE SOP CLASSES SUPPORTED BY STORAGE-SCU**

SOP Class Name	SOP Class UID
12-lead ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1
Ambulatory ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.3
Basic Text SR	1.2.840.10008.5.1.4.1.1.88.11
Basic Voice Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.1
Blending Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.4
Cardiac Electrophysiology Waveform Storage	1.2.840.10008.5.1.4.1.1.9.3.1
Chest CAD SR	1.2.840.10008.5.1.4.1.1.88.65
Color Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.2
Comprehensive SR	1.2.840.10008.5.1.4.1.1.88.33
Deformable Spatial Registration Storage	1.2.840.10008.5.1.4.1.1.66.3
Encapsulated CDA Storage	1.2.840.10008.5.1.4.1.1.104.2
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1
Enhanced SR	1.2.840.10008.5.1.4.1.1.88.22
General ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.2
Grayscale Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.1
Hemodynamic Waveform Storage	1.2.840.10008.5.1.4.1.1.9.2.1
Key Object Selection Document	1.2.840.10008.5.1.4.1.1.88.59
Mammography CAD SR	1.2.840.10008.5.1.4.1.1.88.50
MR Spectroscopy Storage	1.2.840.10008.5.1.4.1.1.4.2
Procedure Log Storage	1.2.840.10008.5.1.4.1.1.88.40
Pseudo-Color Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.3
Raw Data Storage	1.2.840.10008.5.1.4.1.1.66
Real World Value Mapping Storage	1.2.840.10008.5.1.4.1.1.67
RT Beams Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.4
RT Brachy Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.6
RT Dose Storage	1.2.840.10008.5.1.4.1.1.481.2
RT Ion Beams Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.9
RT Ion Plan Storage	1.2.840.10008.5.1.4.1.1.481.8
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5
RT Structure Set Storage	1.2.840.10008.5.1.4.1.1.481.3
RT Treatment Summary Record Storage	1.2.840.10008.5.1.4.1.1.481.7
Segmentation Storage	1.2.840.10008.5.1.4.1.1.66.4

Spatial Fiducials Storage	1.2.840.10008.5.1.4.1.1.66.2
Spatial Registration Storage	1.2.840.10008.5.1.4.1.1.66.1
Standalone Curve Storage (Retired)	1.2.840.10008.5.1.4.1.1.9
Standalone Modality LUT Storage (Retired)	1.2.840.10008.5.1.4.1.1.10
Standalone Overlay Storage (Retired)	1.2.840.10008.5.1.4.1.1.8
Standalone PET Curve Storage (Retired)	1.2.840.10008.5.1.4.1.1.129
Standalone VOI LUT Storage (Retired)	1.2.840.10008.5.1.4.1.1.11
Stereometric Relationship Storage	1.2.840.10008.5.1.4.1.1.77.1.5.3
Stored Print Storage SOP Class (Retired)	1.2.840.10008.5.1.1.27
X-Ray Radiation Dose SR	1.2.840.10008.5.1.4.1.1.88.67

4.2.5.2 Association Policies

4.2.5.2.1 General

STORAGE-SCU initiates but never accepts associations.

**Table 4.2-32
MAXIMUM PDU SIZE RECEIVED AS A SCP FOR STORAGE-SCU**

Maximum PDU size received	114kB (approx.)
---------------------------	-----------------

4.2.5.2.2 Number of Associations

**Table 4.2-33
NUMBER OF ASSOCIATIONS AS A SCP FOR STORAGE-SCU**

Maximum number of simultaneous associations	Unlimited
---	-----------

4.2.5.2.3 Asynchronous Nature

STORAGE-SCU will only allow a single outstanding operation on an Association. Therefore, STORAGE-SCU will not perform asynchronous operations window negotiation.

4.2.5.2.4 Implementation Identifying Information

**Table 4.2-34
DICOM IMPLEMENTATION CLASS AND VERSION FOR STORAGE-SCU**

Implementation Class UID	1.3.6.1.4.1.25403.1.1.1
Implementation Version Name	Dicom 0.1

4.2.5.3 Association Initiation Policy

STORAGE-SCU attempts to initiate a new association for each study, or group of studies, selected by the user (e.g. one association per user-initiated send operation). When initiated by MOVE-SCP, one association is initiated per move request. For move requests initiated from the web interface, one association is used per study. For move requests initiated by auto-routing rules, multiple associations may be used to send a single study.

4.2.5.3.1 Activity – Send Storage Request

4.2.5.3.1.1 Description and Sequencing of Activities

For each instance selected from the user interface to be transferred, a single attempt will be made to transmit it to the selected remote AE. If the send fails, for whatever reason, no retry will be performed, and an attempt will be made to send the next instance.

4.2.5.3.1.2 Proposed Presentation Contexts

**Table 4.2-35
PROPOSED PRESENTATION CONTEXTS FOR STORAGE-SCU**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name	UID		
See Table 4.2-30 (Image SOP Classes)	See Table 4.2-30 (Image SOP Classes)	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCU	None
		RLE Lossless	1.2.840.10008.1.2.5	SCU	None
		JPEG Baseline (Process 1):	1.2.840.10008.1.2.4.50	SCU	None
		JPEG Extended (Process 2 & 4)	1.2.840.10008.1.2.4.51	SCU	None
		JPEG Lossless, Non-Hierarchical (Process 14)	1.2.840.10008.1.2.4.57	SCU	None
		JPEG 2000 Image Compression (Lossless Only)	1.2.840.10008.1.2.4.90	SCU	None
		JPEG 2000 Image Compression	1.2.840.10008.1.2.4.91	SCU	None
See Table 4.2-31 (Non-Image SOP Classes)	See Table 4.2-31 (Non-Image SOP Classes)	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCU	None

STORAGE-SCU will propose Presentation Contexts only for the SOP Class(es) of the instances that are to be transferred.

For each SOP Class being transferred, STORAGE-SCU will propose two Presentation Contexts. One containing the transfer syntax that the SOP class is encoded on for the server, and a second Presentation Context using the uncompressed transfer syntaxes. When the server cannot decompress an image because it does not have the appropriate codec(s), the second Presentation Context (uncompressed syntaxes) is not proposed.

4.2.5.3.1.2.1 Extended Negotiation

No extended negotiation is performed.

4.2.5.3.1.3 SOP Specific Conformance

4.2.5.3.1.3.1 SOP Specific Conformance to Storage SOP Classes

STORAGE-SCU provides standard conformance to the Storage Service Class.

4.2.5.3.1.3.2 Presentation Context Acceptance Criterion

STORAGE-SCU does not accept associations.

4.2.5.3.1.3.3 Transfer Syntax Selection Policies

STORAGE-SCU prefers to send each instance using its current Transfer Syntax, and will find the first Presentation Context where the Transfer Syntax was accepted. In the case where the Transfer Syntax was not accepted by the remote STORAGE-SCP, STORAGE-SCU will check for the acceptance of Explicit VR Little Endian and Implicit VR Little Endian, in that order. If STORAGE-SCU cannot change the Transfer Syntax, the sub-operation will fail and it will not store the instance.

Response Status

STORAGE-SCU will behave as described in the Table below in response to the status returned in the C-STORE response command message.

**Table 4.2-36
RESPONSE STATUS FOR STORAGE-SCU AND RECEIVE STORAGE REQUEST**

Service Status	Further Meaning	Status Codes	Behavior
Failure	Refused: Out of Resources	A7xx	Status logged.
	Data Set does not match SOP Class	A9xx	Status logged.
	Cannot understand	Cxxx	Status logged.
Warning	Coercion of Data Elements	B000	Status logged.
	Data Set does not match SOP Class	B007	Status logged.
	Elements Discarded	B006	Status logged.
Success		0000	Ignored

4.2.5.4 Association Acceptance Policy

STORAGE-SCU does not accept associations.

4.3 NETWORK INTERFACES

4.3.1 Physical Network Interface

The application is indifferent to the physical medium over which TCP/IP executes; which is dependent on the underlying operating system and hardware.

4.3.2 Additional Protocols

When host names rather than IP addresses are used in the configuration properties to specify presentation addresses for remote AEs, the application is dependent on the name resolution mechanism of the underlying operating system.

The ImageServer supports a limited subset of DICOM Web Access to DICOM Objects (WADO). The following MIME types are supported:

**Table 4.3-1
WADO MIME TYPES**

application/dicom	The entire DICOM object is returned.
application/clearcanvas	The raw pixel data as encoded in the DICOM object on the server is returned. This may be compressed or uncompressed pixel data.

The following table contains the WADO parameters that are supported by the ImageServer WADO server.

**Table 4.3-2
WADO PARAMETERS SUPPORTED**

requestType	Must be WADO
studyUID	The Study Instance UID of the requested image.
seriesUid	The Series Instance UID of the requested image.
objectUid	The SOP Instance UID of the requested image.
frameNumber	Only supported for application/clearcanvas MIME type.
contentType	The MIME types as specified in Table 4.3-1.

The URL for accessing the WADO server is structured as follows:

`http://<host>:1000/wado/<serverAE>?requestType=WADO`

By default, the WADO server listens on port 1000. Note also that the AE title of the Server Partition being connected to must be included in the URL.

4.3.3 IPv4 and IPv6 Support

By default, this product supports IPv4. When configured, it will also support IPv6. It does not utilize any of the optional configuration identification or security features of IPv6.

4.4 CONFIGURATION

All configuration is performed through the use of configuration files stored in pre-defined locations that are specific to the underlying operating system.

4.4.1 AE Title/Presentation Address Mapping

The Calling AE Title of the local application is configurable in the configuration file. The mapping of the logical name by which remote AEs are described in the user interface to Called AE Titles as well as presentation address (hostname or IP address and port number) is configurable in the configuration file.

4.4.2 Parameters

**Table 4.4-1
CONFIGURATION PARAMETERS TABLE**

Parameter	Configurable	Default Value
General Parameters		
PDU size	No	N/A
Time-out waiting for acceptance or rejection Response to an Association Open Request. (Application Level timeout)	Yes	10 seconds
General DIMSE level time-out values	No	30 seconds
Time-out waiting for response to TCP/IP connect() request. (Low-level timeout)	Yes	10 seconds
Time-out waiting for acceptance of a TCP/IP message over the network. (Low-level timeout)	Yes	30 seconds
Time-out for waiting for data between TCP/IP packets. (Low-level timeout)	Yes	30 seconds
Send Buffer Size (TCP/IP socket)	Yes	118341 bytes
Receive Buffer Size (TCP/IP socket)	Yes	118341 bytes
Any changes to default TCP/IP settings, such as configurable stack parameters.	No	None
AE Specific Parameters (all AEs)		
Size constraint in maximum object size	No	None
Maximum PDU size the AE can receive	Yes	116794 bytes
Maximum PDU size the AE can send	Yes	116794 bytes
AE specific DIMSE level time-out values	No	30 seconds
Number of simultaneous Associations by Service and/or SOP Class	N/A	N/A
SOP Class support	No	All supported SOP Classes always proposed and accepted
Transfer Syntax support	No	All supported Transfer Syntaxes always proposed and accepted
Other parameters that are configurable	No	None

5. MEDIA INTERCHANGE

The ImageServer does not support any form of Media Interchange.

6. SUPPORT OF CHARACTER SETS

6.1 OVERVIEW

The application supports all extended character sets defined in the DICOM 2008 standard, including single-byte and multi-byte character sets as well as code extension techniques using ISO 2022 escapes.

Support extends to correctly decoding and displaying the correct symbol for all names and strings received over the network, and in the local database.

No specific support for sorting of strings other than in the default character set is provided in the Web GUI. All strings are converted to Unicode before display.

6.2 CHARACTER SETS

In addition to the default character repertoire, the Defined Terms for Specific Character Set in Table 6.2-1 are supported:

Table 6.2-1
SUPPORTED SPECIFIC CHARACTER SET DEFINED TERMS

Character Set Description	Defined Term
Latin alphabet No. 1	ISO_IR 100
Latin alphabet No. 2	ISO_IR 101
Latin alphabet No. 3	ISO_IR 109
Latin alphabet No. 4	ISO_IR 110
Cyrillic	ISO_IR 144
Arabic	ISO_IR 127
Greek	ISO_IR 126
Hebrew	ISO_IR 138
Latin alphabet No. 5	ISO_IR 148
Japanese	ISO_IR 13
Thai	ISO_IR 166
Unicode in UTF-8	ISO-IR 192
Default repertoire	ISO 2022 IR 6
Latin alphabet No. 1	ISO 2022 IR 100
Latin alphabet No. 2	ISO 2022 IR 101
Latin alphabet No. 3	ISO 2022 IR 109
Latin alphabet No. 4	ISO 2022 IR 110
Cyrillic	ISO 2022 IR 144
Arabic	ISO 2022 IR 127
Greek	ISO 2022 IR 126
Hebrew	ISO 2022 IR 138
Latin alphabet No. 5	ISO 2022 IR 148
Thai	ISO 2022 IR 166
Japanese	ISO 2022 IR 13

Japanese	ISO 2022 IR 87
Japanese	ISO 2022 IR 159
Korean	ISO 2022 IR 149
Chinese (Simplified) Extended	GB18030

6.3 CHARACTER SET CONFIGURATION

Whether or not characters are displayed correctly depends on the presence of font support in the underlying operating system. Typically it may be necessary for the user to add one of the “all Unicode” fonts to their system configuration in order to correctly display characters that would not typically be used in the default locale.

7. SECURITY

7.1 SECURITY PROFILES

None supported.

7.2 ASSOCIATION LEVEL SECURITY

None supported.

The system can be configured to only allow configured AEs to open an association.

7.3 APPLICATION LEVEL SECURITY

None supported.

8. ANNEXES

8.1 IOD CONTENTS

8.1.1 Coerced/Modified fields

The following fields may be coerced, based on user input.

**Table 8.1-3
SIGNIFICANT ELEMENTS IN EXPORTED COMPOSITE SOP INSTANCES**

Module	Attribute Name	Tag ID	Value
Patient	Patient Name	(0010,0010)	STORAGE-SCP checks consistency of the tag with other images in the study and will coerce the value when specified through the application Web GUI.
	Patient ID	(0010,0020)	STORAGE-SCP checks consistency of the tag with other images in the study and will coerce the value when specified through the application Web GUI.
	Patient's Birth Date	(0010,0030)	STORAGE-SCP checks consistency of the tag with other images in the study and will coerce the value when specified through the application Web GUI.
	Patient's Sex	(0010,0040)	STORAGE-SCP checks consistency of the tag with other images in the study and will coerce the value when specified through the application Web GUI.
	Issuer of Patient ID	(0010,0021)	STORAGE-SCP checks consistency of the tag with other images in the study and will coerce the value when specified through the application Web GUI.
General Study	Accession Number	(0008,0050)	STORAGE-SCP checks consistency of the tag with other images in the study and will coerce the value when specified through the application Web GUI.

8.2 DATA DICTIONARY OF PRIVATE ATTRIBUTES

No private attributes are defined.

8.3 CODED TERMINOLOGY AND TEMPLATES

Not applicable.

8.4 GRAYSCALE IMAGE CONSISTENCY

Not applicable.

8.5 STANDARD EXTENDED/SPECIALIZED/PRIVATE SOP CLASSES

None

8.6 PRIVATE TRANSFER SYNTAXES

None.